

# BUGTAG™

## Sterile Insect Technique Irradiation Indicator


*It's simple, and can tell you more!*

The BugTag™ SIT series indicators are radiation process verification indicator devices used in Sterile Insect Technique (SIT) programs.



**The BugTag™ SIT Irradiation Indicator** - an innovative irradiation indicator that can tell you more than just "yes" or "no". It can also indicate that the product received at least the minimum dose required when using the color reference system.

**To use:** For a "yes/no" answer, simply confirm that the sensitive dot has turned blue. To verify that the minimum dose has been achieved, simply confirm that the blue dot is the same as or darker than the minimum reference colour. It's that easy!!



### Features and Benefits

- ◆ **Utilizes** a unique colour change material
- ◆ **Quick** and easy to interpret - a simple yes/no process verification indicator
- ◆ **Semi-quantitative** indicator of radiation dose utilizing the minimum reference colour on the label

### Ordering Information

Catalog #:	Use:
BT70	70 Gy
BT90	90 Gy
BT125	125 Gy
BT150	150 Gy
BT200	200 Gy
BT250	250 Gy

**Indicators can also be customized to meet the design and minimum dose needs of your facility**



**To order, contact RadTag Technologies:**

Toll Free: 1-866-RADTAGS  
Phone: (780) 433-4331  
Fax: (780) 432-3726  
E-mail: sales@radtagtech.com

2403 - 96 Street NW  
Edmonton, Alberta  
Canada T6N 0A7  
Web Site: www.radtagtech.com